New Generation Endo Stapling System

Prorduct Name	Product	Description	Color	Trocar diameter
UNIVERSAL Endo Linear Cutter				
	CEAB-30	Endo Linear Cutter (Standard Handle)		
	CEAB-45	Endo Linear Cutter (Short Handle)		
	CEAB-60	Endo Linear Cutter (Extra Long Handle)		
Articulating Reloading units	CADD-30EVS	Endo Articulating reloading unit 30 mm for Vascular/Medium tissue with compression range from 0.88 to 1.5mm	Brown	12mm
	CADD-30ENTS	Endo Articulating reloading unit 30 mm for Medium/Thick		12mm
	CADD-30EN13	tissue with compression range from 1.5 to 2.25mm	Purple	
	CADD-30EPTS	Endo Articulating reloading unit 30 mm for Extra Thick tissue with compression range from 2.25 to 3 mm	• Black	15mm
		Endo Articulating reloading unit 45 mm for Vascular/ Medium		12mm
	CADD-45EVS	tissue with compression range from 0.88 to 1.5mm	Brown	
	CADD-45ENTS	Endo Articulating reloading unit 45 mm for Medium/Thick tissue with compression range from 1.5 to 2.25mm	Durolo	12mm
			Purple	
	CADD-45EPTS	Endo Articulating reloading unit 45 mm for Extra Thick tissue with compression range from 2.25 to 3 mm	Black	15mm
	CADD-60EVS	Endo Articulating reloading unit 60 mm for Vascular/Medium		12mm
		tissue with compression range from 0.88 to 1.5mm	Brown	
	CADD-60ENTS	Endo Articulating reloading unit 60 mm for Medium/Thick tissue with compression range from 1.5 to 2.25mm	Purple	12mm
		Endo Articulating reloading unit 60 mm for Extra Thick		
	CADD-60EPTS	tissue with compression range from 2.25 to 3 mm	Black	15mm





New Generation Endo Stapling System







TECHNOLOGY P3H

- Lateral diffusion of the fluids
- Stepped cartridge surface
- Better tissue compression

The P₃H technology guarantees better stapling delivering better hemostasis and safety


BENEFITS

- Enhanced anvil delivers better compression & manipulation of the tissue
- Stapling with P₃H technology at the areas of greater tissue tension
- Hemostasis & safety

FEATURES

Stepped cartridge Drop-shaped staple groove Progressive staple pushing technology anvil





TECNOIMAGEN

www.tecnoimagen.com.ar | 0810 333 8273