



**Model: 600ROH**

PRESENTATION	The product is packed sterile in box packaging along with separate sterile disposable delivery system to implant the lens in the eye
MATERIAL	Hydrophilic Acrylic CQ UV
OPTIC TYPE	Single Piece, Biconvex, 360° Square Edge
HAPTIC TYPE	Optimized "C"
OPTIC SIZE	6.00 mm
OVERALL SIZE	12.50 mm
ANGULATION	5°
ACD	4.96
REFRACTIVE INDEX	1.462
ESTIMATED A-CONSTANT	118.0
DIOPTER RANGE	-10.0 D to +30.0 D (With 0.5 D steps from +1.0 D to +30.0 D)
IMPLANTATION SITE	Capsular Bag
STERILIZATION	Steam
SHELF LIFE	4 years from date of manufacture
INDICATIONS	600ROH is indicated for implantation in the Capsular Bag as replacement of Natural Human Crystalline Cataractous Lens to restore the refractive functions required in vision process

\* To get best post-op. outcomes, we recommend surgeons to use their optimized A-Constants.

**TECNOIMAGEN**  
www.tecnoimagen.com.ar | 0810 333 8273