

CE
2460

EYECRYL™ PLUS

ACRYLIC FOLDABLE IOL



biotech

VISIONCARE

EYECRYL™ PLUS

ACRYLIC FOLDABLE IOL

FEATURES


- 6 mm optic for better posterior view
- 5° angulation ensures better contact with posterior capsule
- YAG LASER compatible
- UV absorbing material
- Can be injected through 2.2 mm incision



SPECIFICATIONS

MATERIAL	Hydrophilic Acrylic CQ UV									
OPTIC TYPE	Single Piece, 360° Square Edge									
OPTIC SIZE	6.00 mm									
OVERALL SIZE	12.50 mm									
ANGULATION	5°									
ACD	4.96									
REFRACTIVE INDEX	1.462									
ESTIMATED A-CONSTANT	118.0									
RECOMMENDED OPTICAL A-CONSTANTS	<table border="1"> <tr> <td>SRK - T 118.4</td> <td>SRK - II 118.5</td> <td>Holl 1 Const SF : 1.42</td> </tr> <tr> <td colspan="2">HOFFER Q ACD : 5.19</td> <td>HAIGIS α0:0.94, α1:0.40, α2:0.10</td> </tr> <tr> <td>Barrett 1.52</td> <td>Hill-RBF 118.30</td> <td>Olsen ACD 4.20</td> </tr> </table>	SRK - T 118.4	SRK - II 118.5	Holl 1 Const SF : 1.42	HOFFER Q ACD : 5.19		HAIGIS α0:0.94, α1:0.40, α2:0.10	Barrett 1.52	Hill-RBF 118.30	Olsen ACD 4.20
	SRK - T 118.4	SRK - II 118.5	Holl 1 Const SF : 1.42							
	HOFFER Q ACD : 5.19		HAIGIS α0:0.94, α1:0.40, α2:0.10							
Barrett 1.52	Hill-RBF 118.30	Olsen ACD 4.20								
DIOPTRER RANGE	-10.0 to +30.0 D (with 0.5 D steps from +1.0D to +30.0D)									
IMPLANTATION SITE	Capsular Bag									
STERILIZATION	Steam									
SHELF LIFE	4 years from date of manufacture									

Model: 600



TECNOIMAGEN
www.tecnoimagen.com.ar | 0810 333 8273