



Enquire

# SAF 6125SQ



Optic Diameter	6.00 mm	
Overall Diameter	12.50 mm	
Optic Design	Aberration Neutral, Negative Aspheric Biconvex	
Haptic Design	"C" design, 5° Angulation	
PCO Protection	360° perfect Square edge around the optic edge	
Nominal	A	118.0
Haigis	a0	0.568
Haigis	a1	0.400
Haigis	a2	0.100
HofferQ	pACD	4.80
Holl.	1 Const sf	1.02
SRK/T	A	117.88
SRK II	A	118.0
Refractive Index	1.4570	

## Features

- ▶ Hydrophilic Acrylic, 26% Water content with UV filter
- ▶ UA-A (Violet light and dark blue light) filter and allowing healthy blue spectrum light in
- ▶ Power Available :- 10.0D to +40.0D
- ▶ Better visual acuity and Enhanced contrast even at night mode.(+10.0D to +30.0D in 0.5D steps -10.0D to +9.0D and +31.0 to 40.0D in 1.0D steps