




NEW PANTHER SMART STAPLER



- Smart reload identification
- Smart tissue recognition
- Smart firing control



Panther[®]
Healthcare

ELECTRONICALLY ASSISTED

Chip incorporated in the reloads for the most precise firing

Innovative 45 degrees Non-Stop articulation system

SMART HANDLE

Universal handle ready to use with the most extended portfolio of reloads



One touch firing and opening system

Universal handle with better design and ergonomics for easy 360 degrees knobbing



P3H

TECHNOLOGY P3H

- Lateral diffusion of the fluids
- Stepped cartridge surface
- Better tissue compression

The P3H technology guarantees better stapling delivering better hemostasis and safety

BENEFITS

- Enhanced anvil delivers better compression & manipulation of the tissue
- Stapling with P3H technology at the areas of greater tissue tension
- Hemostasis & safety



Purple Black



Brown



Curved Tip Reloads for the most challenging procedures

FEATURES

Stepped cartridge



Drop-shaped staple groove



Progressive staple pushing technology



Enhanced anvil



Stronger clamping force



Four staples at the end of the staple line


















Interchangeable cartridge



Coating



NEW GENERATION ENDO STAPLING SYSTEM

Product Name	Product	Description	Color	
Smart Stapler	PEAS	Powered Endoscopic Cutter Stapler with Non-Stop articulation system (Short Handle)		
	PEAM	Powered Endoscopic Cutter Stapler with Non-Stop articulation system (Standard Handle)		
	PEAL	Powered Endoscopic Cutter Stapler with Non-Stop articulation system (Extra long Handle)		
Articulating Reloading units	CADD-30EVS	Endo Articulating reloading unit 30 mm for Vascular/Medium tissue with compression range from 0.88 to 1.5mm	 Brown	12mm
	CADD-30ENTS	Endo Articulating reloading unit 30 mm for Medium/Thick tissue with compression range from 1.5 to 2.25mm	 Purple	12mm
	CADD-45EVS	Endo Articulating reloading unit 45 mm for Vascular/ Medium tissue with compression range from 0.88 to 1.5mm	 Brown	12mm
	CADD-45ENTS	Endo Articulating reloading unit 45 mm for Medium/Thick tissue with compression range from 1.5 to 2.25mm	 Purple	12mm
	CADD-45EPTS	Endo Articulating reloading unit 45 mm for Extra Thick tissue with compression range from 2.25 to 3 mm	 Black	15mm
	CADD-60EVS	Endo Articulating reloading unit 60 mm for Vascular/Medium tissue with compression range from 0.88 to 1.5mm	 Brown	12mm
	CADD-60ENTS	Endo Articulating reloading unit 60 mm for Medium/Thick tissue with compression range from 1.5 to 2.25mm	 Purple	12mm
	CADD-60EPTS	Endo Articulating reloading unit 60 mm for Extra Thick tissue with compression range from 2.25 to 3 mm	 Black	15mm
Articulated Reloads for endo linear cutter stapler	CADE-30EVS	Endo Articulating reloading unit 30 mm with Curve Tip and variable-sized staple height 2.0 mm, 2.5 mm, 3.0 mm	 Brown	12mm
	CADE-30ENTS	Endo Articulating reloading unit 30 mm with Curve Tip and variable-sized staple height 3.0 mm, 3.5 mm, 4.0 mm	 Purple	12mm
	CADE-45EVS	Endo Articulating reloading unit 45 mm with Curve Tip and variable-sized staple height 2.0 mm, 2.5 mm, 3.0 mm	 Brown	12mm
	CADE-45ENTS	Endo Articulating reloading unit 45 mm with Curve Tip and variable-sized staple height 3.0 mm, 3.5 mm, 4.0 mm	 Purple	12mm
	CADE-60EVS	Endo Articulating reloading unit 60 mm with Curve Tip and variable-sized staple height 2.0 mm, 2.5 mm, 3.0 mm	 Brown	12mm
	CADE-60ENTS	Endo Articulating reloading unit 60 mm with Curve Tip and variable-sized staple height 3.0 mm, 3.5 mm, 4.0 mm	 Purple	12mm
	CADE-60EPTS	Endo Articulating reloading unit 60 mm with Curve Tip and variable-sized staple height 4.0 mm, 4.5 mm, 5.0 mm	 Black	15mm

TECNOIMAGEN

www.tecnoimagen.com.ar | 0810 333 8273