

|                               | SK-5000A                                                                                 | SK-5000B                                                                                              |
|-------------------------------|------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------|
| Measurement mode              | Fully automatic/Automatic/Manual                                                         |                                                                                                       |
| Corneal thickness measurement | Not support                                                                              | Support<br>Measurement range-CCT:400-800um<br>Measurement increment-CCT:10um<br>Light Source:Blue LED |
| Display                       | 8-inch 1024*600 color touch screen<br>support up and down 180° , left and right 60° flip |                                                                                                       |
| Fixation light                | internal fixation light, support                                                         |                                                                                                       |
| Measurement accuracy          | 1mmHg(0.133Kpa)                                                                          |                                                                                                       |
| Measurement range             | 0-60mmHg (0-8Kpa)<br>Include automatic switching between30/60mmHg                        |                                                                                                       |
| Printer                       | Built-in thermal printer, automatic paper cutting                                        |                                                                                                       |
| Data output/input             | Input:USB; Output:LAN/RS232C                                                             |                                                                                                       |
| Size and weight               | 500*300*500; <25kg                                                                       |                                                                                                       |

**TECNOIMAGEN**  
[www.tecnoimagen.com.ar](http://www.tecnoimagen.com.ar) | 11-4582-2222

QUICK SPEED  
 3D AUTO-TRACKING  
 AUTO MEASUREMENT  
 SOFT AIR PUFF, SAFTY  
 PACHYMETER



**NON-CONTACT  
 TONOMETER**  
 SK-5000A/B





Electric chin rest




SK-5000B supports corneal thickness measurement  
More functions



Auto alignment  
Auto focus




Fully automatic mode  
One touch  
You can get OS&OD data  
It also supports semi-automatic/manual mode



Corneal thickness measurement  
Auto soft air puff



Support air puff detection after turn on  
Support safety limit setting  
Safer equipment use and operation



Auto switch OS&OD



Print report  
Auto paper cutting



8 inch color capacitive touch screen  
Can be flip up and down 180 °  
Flip right and left 60°  
Operation and observation are more convenient



**ONE CLICK**

**FULLY AUTOMATIC**

