




# NEW PANTHER SMART STAPLER



Smart reload identification  
Smart tissue recognition  
Smart firing control



**Panther**<sup>®</sup>  
Healthcare

## ELECTRONICALLY ASSISTED

Chip incorporated in the reloads for the most precise firing

Innovative 45 degrees Non-Stop articulation system

## SMART HANDLE

Universal handle ready to use with the most extended portfolio of reloads



One touch firing and opening system

Universal handle with better design and ergonomics for easy 360 degrees knobbing



# P3H

## TECHNOLOGY P3H

- Lateral diffusion of the fluids
- Stepped cartridge surface
- Better tissue compression

The P3H technology guarantees better stapling delivering better hemostasis and safety

## BENEFITS

- Enhanced anvil delivers better compression & manipulation of the tissue
- Stapling with P3H technology at the areas of greater tissue tension
- Hemostasis & safety



Purple Black



Brown



Curved Tip Reloads for the most challenging procedures

## FEATURES

Stepped cartridge



Drop-shaped staple groove



Progressive staple pushing technology



Enhanced anvil



Stronger clamping force



Four staples at the end of the staple line



Interchangeable cartridge



Coating



**TECNOIMAGEN**

[www.tecnoimagen.com.ar](http://www.tecnoimagen.com.ar) | 11-4582-2222