

mindray

HyLED Q2M
LED Examination Light

Brightness follows your motion



Improving ergonomics

- Handle rotation enables illumination adjustment(optional) with precision control.
- Adjust the surgical light's on/off, illumination intensity via intuitive touch buttons.
- LED service life up to 80,000 hours with 5 years warranty.



Four lockable wheels



Long-lasting backup battery



Handle for adjusting illumination (optional)

Infection control

- The light head's seamless surface minimizes dirt accumulation and cross-contamination risks.
- Validated for cleaning with various disinfectants.
- Detachable sterilizable handle (optional).

Technical Specifications*

Central illuminance (at 1m distance) **	70,000 lux
Light field diameter (at 1m distance)	170mm
Light field (D50/D10)	> 50%
Depth of illumination (20%)	1,500mm
Depth of illumination (60%)	1,100mm
Color temperature	4,350K
Color rendering index (Ra)	99
Color rendering index (R9)	97
Power supply of all light sources	25W

* All values measured according to IEC 60601-2-41.

** Due to manufacturing and measuring tolerances, illuminance has a tolerance of +/- 10%.