



SINOANGIO-1200

High-pressure Injector



Sino Medical-Device Technology Co., Ltd.

📍 6th Floor, Building 15, No.1008.Songbai Road, Nanshan District, Shenzhen, P.R.China Zip: 518055

☎️ +86-755-26755675 +86-755-26755692 Hotline: +86-4008860334

SinoAngio-1200

High-pressure Injector



Feature and Benefits

Safety function	<ul style="list-style-type: none">▪ Unique "Direct Pressure Sensor" detects real-time injection pressure , stop injection automatically if the pressure is abnormal, very helpful to reduce the swelling caused by vascular rupture▪ Stop Injection after alarm if Injection blockage occur▪ "Air exhaustion confirm" button prevent air embolism▪ "Trial injection" button avoid fault operation▪ Soft and hard "stop" button▪ Sound and optic alarm▪ Design of anti-leakage liquid in the head to prevent damage caused by the liquid entering inside of Injector head Real-time display of injection pressure curve during injection
Components quality	<ul style="list-style-type: none">▪ Best quality motor and key parts , long-term working robust.▪ Up to date powerful 32X processor , real-time control
Injection functions	<ul style="list-style-type: none">▪ Flow rate : 0.1-50.0ml/s , increments 0.1ml/s▪ Pressure limits : 100 -1200 Psi▪ Scan delay : 0-600 s▪ injection delay : 0-600 s▪ Max . 8 phases support different injection combination▪ Max . 120 injection programs can be stored , recalled and modified▪ 150ml large volume syringe meets most application needs
Operability	<ul style="list-style-type: none">▪ Local & remote dual monitors , easy to operate▪ 12.1 " true-color touch screen LCD , convenient graphic user interface▪ Extensive in-and output modes can be fitted to various imaging devices▪ Mobile cart convenient for relocation▪ Shine LED indication , show injecting or filling status visually

Specifications

Model	SinoAngio-1200
Syringes	150ML , single syringe
Power	a.c. 220V 50/60 Hz
Weight	67KG
Environment	Temperature:+5°C~+40°C Humidity: ≤ 93%
Dimensions	Main machine:1300mm(Height)×530mm(Width)×530mm(Depth)

Configuration

Main machine , documents package , wire package

Manufacturer : Sino Medical-Device Technology Co., Ltd.

Add : 6th Floor , Building 15 , No.1008 , Songbai Road , Nanshan District , Shenzhen City.P.R.China.
Tel : (86)755-26755692 Fax: (86)755-26755687 Zip:518052
Website : <http://www.sinomdt.com> Email : information@sinomdt.com